

Si-Shift

High-performant non-oriented electrical steels with a silicon content beyond today's limits:
new materials for an electrified future.



Deliverable D2.1 – Mössbauer Spectroscopy Protocol

Author: Elżbieta Jartych, Lublin University of Technology, Poland

Si-Shift - Project: 101112518 —RFCS-02-2022-RPJ

Start date: 2023-09-01

Duration: 42 months

Project coordinator: Ghent University

Coordinator contact: Stefaan Cottenier, Stefaan.cottenier@ugent.be



Funded by
the European Union

Document information

Project Acronym	Si-Shift
Project No.	101112518
Project Call	RFCS-02-2022-RPJ - Steel Research Projects
Title	Mössbauer spectroscopy protocol
Deliverable no.	D2.1
Draft delivery date	30.05.2024 (deadline extended to 31.08.2024)
Dissemination level	PU
Deliverable type	Report
Authors	Elżbieta Jartych

Document history:

Version	Date of issue	Content and changes	Edited by
1.0	30.05.2024	First draft	Elżbieta Jartych
1.1	16.07.2024	Additions: figures with exemplary CEMS spectra	Elżbieta Jartych
1.2	26.08.2024	Additions: description of the software	Elżbieta Jartych



Funded by
the European Union

Abstract

Deliverable abstract

The Mössbauer spectroscopy protocol will deliver all the information necessary to perform an analyze Mössbauer experiments on Fe-Si samples, needed to study an order-disorder transition in Fe-Si system.

The document describes the theoretical background of Mössbauer spectroscopy and details of the experimental setup. It defines the elements and materials necessary for the measurements. The main part of the protocol includes the description of preparing the spectrometer, calibration of the spectrometer transducer speed, software, sample measurement and analysis.

As an illustration, the measurement and analysis of one sample is shown.



Funded by
the European Union

Contents

1	Introduction	5
2	Theoretical background of Conversion Electron Mössbauer Spectroscopy CEMS	5
3	Construction of a spectrometer	8
4	Conversion electron counters	9
5	Materials necessary for CEMS measurements	12
6	CEMS procedure	13
	6.1 Preparing the spectrometer	13
	6.2 Calibration of the spectrometer transducer speed	13
	6.3 Sample measurement	14
	6.4 Software	14
7	Results of measurements	18
8	References	20

1 Introduction

As described in the Si-Shift proposal (1.2.2. Pillar 1: Material exploration), Mössbauer spectroscopy will be used as a tool to get information about local atomic arrangement in Fe-Si and Fe-Si-X samples. The method will complement *ab initio* calculations combined with machine learning and more conventional experimental characterization methods (e.g., X-ray diffraction, transmission electron microscopy TEM, etc.).

Mössbauer spectroscopy is a local technique (i.e., providing information at the atomic scale) using gamma radiation. It is especially suitable for iron-containing crystals because the naturally occurring ^{57}Fe isotope has ideal properties for this technique. Mössbauer spectroscopy provides a fingerprint spectrum with contributions from all types of surroundings of iron atoms in the alloy, at the atomic scale (sensitivity of a few atomic neighbour shells). It is possible to numerically deconvolute an experimental Mössbauer spectrum into single components, each reflecting a specific local configuration of atoms.

In the Si-Shift project, all the investigated samples are Fe-Si electrical steel in bulk form. At the stage of creating the Si-Shift proposal, it was established that it is impossible to produce thin sheets of these alloys with a thickness of 40-50 micrometers. Therefore, the Conversion Electron Mössbauer spectroscopy CEMS was established as a suitable technique for studies in the project.

2 Theoretical background of Conversion Electron Mössbauer spectroscopy CEMS

The effect of recoilless emission and absorption of γ radiation by the nuclei in solids discovered in 1957 by a German physicist, Rudolf L. Mössbauer became the basis for the development of a research technique called Mössbauer Spectroscopy (MS). The most popular Mössbauer nuclides are ^{57}Fe , ^{119}Sn , ^{151}Eu , ^{155}Gd , and ^{161}Dy because they exhibit the largest f factor, which means the probability of the Mössbauer effect.

It is the ratio of the recoilless photons that are emitted to the total number of photons, called also the Lamb-Mössbauer factor:

$$f = \exp(-k^2 \langle x^2 \rangle), \quad (1.1)$$

where $k = 2\pi/\lambda$ is a wave vector and $\langle x^2 \rangle$ denotes the square average value of the amplitude of atom oscillations in the direction of the quantum momentum that is emitted. The value of the f factor strongly depends on the temperature and the energy of the transition of the nucleus from an excited to a ground state. Thus, the probability of the Mössbauer effect increases with a decrease in the measurement temperature and is larger for low energy transitions. In the case of isotope ^{57}Fe , the value of $f = 0.75$ at room-temperature for energy $E = 14.4$ keV. The iron is prevalent in nature and materials, and in all of the samples that contain iron, the relative concentration of ^{57}Fe isotope is equal to 2.17 %.

The main idea of the Mössbauer spectroscopy is the motion of the source of γ radiation (Figure 1). The movement of the source relative to the absorber (the investigated sample) causes a change in the energy of γ quanta in accordance with the Doppler effect:

$$\Delta E = (V/c) E_{\gamma}, \quad (1.2)$$

where E_{γ} is the energy of the γ quantum and V denotes the velocity of the source relative to the absorber. For γ radiation with an energy of 14.4 keV, the change of velocity by 1 mm/s causes a modulation of energy by 48 neV. This modulation allows the Mössbauer spectrum to be registered, i.e., the possible shifts and splitting of the absorption lines can be observed. The sensitivity of Mössbauer spectroscopy is one of the highest in experimental physics. The uncertainty of the determination of energy for ^{57}Fe is about 10^{-10} eV. The most important problem is to precisely determine the movement of the source with constant acceleration. The range of the source motion that is selected is dependent on the composition of the material.




Figure 1. Schematic diagram of an experimental setup for Mössbauer spectroscopy in transmission geometry.

Two popular techniques are used to measure Mössbauer spectra. The first one is a transmission geometry in which the sample is placed between the mobile source and the detector (Figure 1). In this case, the absorber should be made of tens mg/cm^2 thick foil. The most effective technique for investigating thin films and surface layers of materials is Conversion Electron Mössbauer Spectroscopy – CEMS in which the sample is placed inside a conversion electrons counter that is filled with a mixture of gases. In both techniques, transmission Mössbauer spectroscopy (TMS) and Conversion Electron Mössbauer Spectroscopy (CEMS), the same recoil-free nuclear resonance is the fundamental principle of the method; however, a different product of the resonance is utilized. Instead of 14.4 keV photons (in the case of TMS based on ^{57}Fe), 7.3 keV conversion electrons together with 5.4 and 6.2 keV Auger electrons are detected in CEMS technique (Figure 2).




Figure 2. De-excitation of ^{57}Fe nucleus after resonance excitation [1].

In comparison with TMS, the CEMS technique does not require the minimum absorber thickness due to the much lower detection limit of CEMS. Moreover, the internal conversion process is much more probable than the emission of γ radiation. The fundamental principle of the CEMS measurement is presented in Figure 3. The penetration depth of 14.4 keV photons is about 40 μm , while that of 7.3 keV electrons reaches up to 350 nm. If the resonance occurs in bulk (Figure 3 a), all electrons are absorbed by a sample material, and they are no longer registered. However, when the resonance occurs in the upper surface layer, some electrons can escape and be detected (Figure 3 b) [2].




Figure 3. The principle of the CEMS measurement [2].

Mössbauer spectroscopy is a technique that can be used to investigate hyperfine interactions. These are interactions of the electromagnetic moments of the nuclear probe (i.e., ^{57}Fe isotope) with the electric and magnetic fields that originate from its electron shells as well as from the nearest neighborhood of the probe. A detailed description of hyperfine interactions can be found in many monographs and scientific articles (e.g., [3-5]). There are three the most important interactions: (i) the monopole electrostatic interaction between the positive charge of the atomic nucleus and the s -type electrons, (ii) the interaction between the quadrupole electric moment of the nucleus and the electric field gradient that is generated by the surrounding electrons and (iii) the interaction between the magnetic moment of the nucleus and the magnetic field at the nucleus site.

In general, all of the interactions mentioned above can occur simultaneously in the investigated material. The possible shifts and splitting of nuclear levels are reflected in the Mössbauer spectrum that is registered. The numerical fitting of the experimental spectrum together with the structural information, for example from the XRD patterns, allows the phase composition of the material to be recognized and its magnetic properties to be determined. The fitting procedure is performed until the best fitting parameter, χ^2 , is achieved. The hyperfine interaction parameters that are determined from the spectra are as follows: IS – the isomer shift relative to α -iron, QS – the quadrupole shift of the magnetically split spectra or the quadrupole splitting of the paramagnetic doublet, B_{hf} – hyperfine magnetic field, Γ – half width at half maximum of the spectral lines. All of the parameters are determined with an accuracy of about 0.01 mms^{-1} , which translates into an accuracy of the determination of the hyperfine interaction energy of the order of 10^{-10} eV .

3 Construction of a spectrometer

The block diagram of the CEMS spectrometer is presented in Figure 4.




Figure 4. The CEMS measurement setup.

The heart of the system is a head with a moving radioactive source. The movement of the source causes a tiny modulation of energy of outgoing gamma quanta due to the Doppler effect. The source of 14.4 keV gamma radiation is ^{57}Co embedded in a cubic and paramagnetic matrix to avoid hyperfine structure of the source being overlaid on a hyperfine structure of the sample for investigation. The matrix usually is Rh or Cr. The source is attached to the end of the core of an electro-mechanic transducer. The coil of the transducer is supplied from a function generator. Various shapes of the signal can be used to drive the transducer, according to the needs of a user; however, a triangular signal leading to the constant acceleration working mode of a spectrometer is used most frequently.

The sample for investigations is located inside the detector of electrons and grounded. High voltage from the power supply is applied between the sample and detection rods playing the role of a positive electrode. Thus, when the gamma quanta hitting the sample causes the emission of internal conversion electrons, they are accelerated towards positively charged detection rods. In this way electric pulses are generated as a consequence of conversion electron emission. Additionally, a very slow gas flow (~ 0.1 l/h) is necessary to ensure the stable work of the detector. The spike pulses coming from the detector are amplified and formed. To cut the low-amplitude electronic noise the discriminator is applied. The latter

transforms spike pulses to TTL pulses that can be counted by a main multichannel analyzer working in a time domain (TMA). There is also another, auxiliary MA working in an amplitude domain (AMA). The amplitude spectrum of pulses obtained from the latter one can be used for discrimination threshold regulation by the user. The TMA is synchronized with the function generator driving the transducer thus the coming TTL pulses can be categorized into channels regarding the time of coming. After a calibration process, a specific velocity/energy can be ascribed to each time channel. In this way, a velocity/energy spectrum of gamma quanta causing the internal conversion phenomenon can be created on a PC.

4 Conversion electron counters

Because the probability of the emission of conversion electrons from the K shell is the largest (77 – 81 %), mainly electrons with low energy, i.e., 7.3 keV are detected. However, the detection of K-electrons is disturbed by an intensive gamma background. In order to suppress the gamma photons, a gas-flow proportional counter can be used in which helium is a counting gas with an admixture of methane as a quencher. The most important requirements for the counter are long-term stability and suppression of non-resonant background (with the minimum self-produced photoelectrons). Further parameters to consider are detector thickness, number of wires, wire diameter, gas mixture, gas-flow rate, and fixing the sample [6]. The scheme of a simple gas-flow proportional detector for CEMS measurements is presented in Figure 5.




Figure 5. A simple gas-flow proportional detector for CEMS [2].

To suppress the gamma background, helium is used as a detection gas with an admixture of 4 – 12 % methane as a quenching gas. In the dependence of the methane concentration, there are various widths of the proportional band, i.e., between 1000 and 1360 V for 8 % of CH_4 and between 1145 – 1540 V for 15 % CH_4 [2]. The wider proportional band, the measurement is more stable, especially in the case of long-lasting measurements. The amount of gas flow between a few and a hundred cm^3 per hour is enough for good measurement.

As an example, a description for two conversion electron counters with different designs is presented as follows.

Counter LEK 2

The actual appearance of the LEK 2 detector is shown in Figure 6.



Figure 6. The LEK 2 conversion electron detector.

In this type of the counter, the sample must have certain dimensions. Area of the sample is $\sim 2 - 4 \text{ cm}^2$ (in the form of sheet $2 \text{ cm} \times 2 \text{ cm}$ or disc with diameter 2 cm). The maximum thickness of the sample is less than $\sim 5 \text{ mm}$. The sample is mounted inside the counter. Two electrodes collect charge, which is directed to the preamplifier.

Counter MCD-14

The actual appearance of the MCD-14 detector is shown in Figures 7 and 8.




Figure 7. The MCD-14 conversion electron detector.




Figure 8. The elements of MCD-14 conversion electron detector and a chamber with four electrodes.

In the MCD-14 detector, the area of the sample is still $\sim 2 - 4 \text{ cm}^2$ (in the form of sheet $2 \text{ cm} \times 2 \text{ cm}$ or disc with diameter 2 cm); however, the maximum thickness of the sample is less than $\sim 3 \text{ mm}$. The sample is mounted inside the counter and there is a possibility to change the distance between the sample and electrodes. Four electrodes collect charge, the efficiency of counting is higher than in the case of LEK 2 counter.

The real view of the mounted MCD-14 detector is presented in Figure 9.




Figure 9. The experimental setup with MCD-14 counter.

5 Materials necessary for CEMS measurements

Necessary for the CEMS measurements is the radioactive source, i.e., nuclide that emits γ quanta with the suitable energy. In the case of the samples containing iron, the proper source is ^{57}Co which decays during electron capture to the ^{57}Fe stable isotope (scheme of the decay is shown in Figure 10).




Figure 10. The scheme of the ^{57}Co radioactive decay.

In measurements, the energy of γ quanta 14.4 keV is used. The activity of the source about 1 GBq is enough to lead the effective measurements for 3 years. The sketch of the source ^{57}Co in Rh matrix is presented in Figure 11. The dimensions given are standard and fit all heads in Mössbauer spectrometers used in laboratories around the world. The minimal activity of the source necessary for measurements is about 1 GBq.




Figure 11. The sketch of the source ^{57}Co with cover and mounting element (M4).

The second important element is the mixture of gases. In the case of composition 90 % helium and 10 % methane, the supply voltage on the electron counter is typically in the range of 1000 – 1300 V. The volume of 50 liters is enough to carry out measurements by 1 year.

6 CEMS procedure

The procedure of CEMS measurements includes three main steps.

6.1 Preparing the spectrometer

To prepare the spectrometer to measurements is needed:

- installation of a reference (calibration) sample - α -Fe foil with a thickness of 25 μm ;
- installation of an internal conversion electron detector in the spectrometer;
- regulation of the optimal gas flow through the detector (90 % He + 10 % CH_4 mixture);
- recording a continuous energy spectrum of secondary electrons resulting from the phenomenon of the internal conversion, observing electrical impulses on an oscilloscope;
- selection/correction of operating parameters based on the observed pulse spectrum (electron energy spectrum):
 - supply voltage of the electron counter (voltage between the sample and the electrodes), typically in the range 1000 – 1300 V;
 - gain of the amplitude of electrical impulses from the counter, usually 40 – 100 times;
 - pulse width 0.5 – 2 μs ,
 - selecting the pulse discrimination threshold - cutting off electronic noise; the discrimination threshold is usually in the range of 50 – 100 mV depending on the set gain.

6.2 Calibration of the spectrometer transducer velocity

In the case of measurements in constant acceleration mode:

- a) selection of the maximum velocity (amplitude of the triangular signal) adequate to the expected spectrum of the tested sample;
- b) starting the spectrometer and recording the calibration spectrum; for a radioactive source with an activity of 1 GBq collecting time is between 2 and 6 h;
- c) automatic calibration – assigning a speed (energy) value to each time channel in which the pulses were recorded (it is performed automatically by the software).

6.3 Sample measurement

Main sample measurement includes:

- a) replacing the calibration sample by the test sample; the sample should be conductive and have a thickness proper to the type of counter;
- b) selecting the appropriate distance between the sample and the detector electrodes (according to the sample thickness);
- c) starting the measurement with the spectrometer operating parameters identical to those of the calibration sample; the measurement lasts from 1 day to 2 weeks depending on the sample type and iron content;
- d) recording the Mössbauer spectrum with the appropriate number of detection channels; the maximum possible number of channels is 2×2048 ; this choice compromises between a sufficiently large counting statistic and a sufficiently small speed step (energy resolution); typically, the spectrum is registered in 2×512 channels.

6.4 Software

The software used, program MOS, was elaborated by Jan Żukrowski on the basis of the papers [7-11]. The numerical program allows fitting single lines, doublets and Zeeman sextets depending on the kind of Mössbauer spectrum. The maximum number of the components is 36.

As a results of fitting, the following parameters are obtained for the suitable component:

- hyperfine interactions parameters – isomer shift, quadrupole splitting (shift), hyperfine magnetic field;
- spectral parameters – amplitude of lines, half width at half maximum, ratios intensity of spectral lines 1 to 3 and 2 to 3, area of spectral lines.

The Lorentzian shape of spectral lines is assumed during the fitting procedure.

The MOS software works in an iterative manner. It allows one to obtain the best set of hyperfine interactions and spectral parameters, i.e., the set giving the lowest χ^2 parameter.

Experimental data can be loaded in two types of files:

- a) an 8-column text file with the number of counts in each channel separated by tab (Figure 12);
- b) a 2-column text file with the first column being velocity and, the second column being the counts number for each channel (Figure 13).

The starting hyperfine parameters for the fitting procedure are gathered at the end of the file as shown in Figure 12. Besides that, just below the experimental data, number of channels, folding point, and maximum velocity must be included.

cech_19_08_2020_fit

X

cech_19_08_2020_fit_2columns

+

Plik

Edytuj

Wyświetl

19.08.2020

cechowanie z g₁ 5.0 matr Rh nowe 512 zlo

512

20622	20678	20684	20791	20757	20844	20771	20788
20856	20811	20946	20613	20835	21054	20979	20930
20883	20980	21033	20633	20914	20876	20666	20867
20896	20647	20971	20775	20694	20945	20723	21041
21034	20784	20825	20807	20788	20951	20775	20735
20944	20691	20737	20874	20917	20842	20902	20716
21100	20652	20730	20862	20802	20913	20716	21002
20558	20808	20760	20811	20483	20588	20796	20806
20799	20843	20842	20745	20882	20810	20913	20806
20919	20805	20677	20597	20793	20954	20796	20638
20716	20837	20928	20660	20938	20802	20753	20750
20654	20656	20863	21072	20842	20681	20800	20882
21067	20789	20721	20635	20436	20658	20724	20828
20834	20780	20820	20810	20766	20824	20629	20687
20751	20632	20689	20638	20580	20485	20413	20670
20562	20672	20300	20152	20119	20329	20070	19647
19419	18903	17988	16863	16082	15635	15701	16315
17145	18162	18723	19381	19744	19928	19941	19958
19966	20100	20626	20320	20678	20441	20531	20464
20691	20439	20715	20584	20667	21005	20805	20490
20702	20901	20737	20566	20561	20328	20612	20814

20240	20028	20431	19972	19874	19677	19405	19228
18897	18293	17535	16555	15503	15314	15821	16835
17875	18644	19368	19576	19807	19999	20232	20181
20431	20190	20611	20531	20324	20492	20649	20616
20419	20581	20881	20583	20568	20843	20611	20907
20685	20570	20838	20803	20722	20818	20874	20579
21043	20935	21019	20889	20844	20844	20692	20647
20614	20701	20931	20540	20709	20779	20847	20850
20688	20837	20998	20578	20752	20674	20867	20819
20869	20743	20838	21000	20787	20849	21040	20725
20822	20858	20972	20849	20928	20873	20689	20781
20482	20891	21042	20967	21052	20831	20683	20611
20871	20799	20889	20727	20803	20795	20604	20766
21138	20660	20617	20552	20721	20675	20746	20815
20719	20566	20778	20921	20899	20750	20706	20813
20864	21230	20985	20817	21011	20873	20973	21105
21077	20967	20941	20781	20798	20832	20873	20890
20888	20892	20741	20685	20752	20714	20966	20654
20804	20904	20632	20648	21069	20674	20872	20775

1024

512.5964

1

11.3640

'Fe'

0

1

512

1|

1

'B1'

20748

0

0

1

'B2'

0.0

1

1

1

'B3'

0.0

1

1

1

'FSB'

0.0

1

1

2

'C'

12774

0

0

1

3

'IS'

-0.10

0

0

1

4

'H'

332.76

0

0

1

5

'QS'

0.00

0

0

1

6

'G1'

0.142

0

0

1

6

'G2'

0.137

0

0

1

6

'G3'

0.126

0

0

1

7

'CL2'

2.0

0

0

1

Figure 12. The input file with experimental data and starting hyperfine parameters in the 8-column text file.

cech_19_08_2020_fit			
cech_19_08_2020_fit_2columns			
Plik	Edytuj	Wyświetl	
19.08.2020 cechowanie z gł 5.0 matr Rh nowe 512 zło			
512			
1024	512.5964	1	11.36400 'primary data number, folding point, velocity code, Vmax
-11.342	20622		
-11.297	20678		
-11.253	20684		
-11.209	20791		
-11.164	20757		
-11.120	20844		
-11.075	20771		
-11.031	20788		
-10.987	20856		
-10.942	20811		
-10.898	20946		
-10.854	20613		
-10.809	20835		
-10.765	21054		
-10.720	20979		
-10.676	20930		
-10.632	20883		
-10.587	20980		
-10.543	21033		
-10.498	20633		
-10.454	20914		
-10.410	20876		
-10.365	20666		
-10.321	20867		
-10.276	20896		
-10.232	20647		
-10.188	20971		
-10.143	20775		
-10.099	20694		
-10.054	20945		
-10.010	20723		
-9.966	21041		
-9.921	21034		
-9.877	20784		
-9.833	20825		
-9.788	20807		
-9.744	20788		

Figure 13. Experimental data in the 2-column text file.

The latest version of the MOS software is MosK program. The advantage of the MosK software is the possibility of pre-processing of experimental data (Figure 14). Among the most useful operations possible in this step are the following:

- folding of the unfolded spectrum, either in automatic or manual mode;
- change of number of channels for the registered spectrum. This step allows one to choose the best channel number, i.e., the number that gives the best signal-to-noise ratio with acceptable energetic resolution;
- adding two spectra in the case of 2-spectra were being registered simultaneously, e.g., by use of 14.4 keV line and the "escape peak";
- including calibration data;
- closer inspection of the spectrum in the "point-by-point" manner.




Figure 14. The front panel of the MosK software.

Loading of the input file is executed by clicking file -> OpenSourceFile (Figure 15).




Figure 15. The page for Source File loading in the MosK ver. 12.2010 software.

After loading the source file, the window with the experimental and calculated pattern is opened (Figure 16). The last step is an iterative spectrum fitting using the internal procedure based on the χ^2 minimization. After this step, the user must check whether the procedure was convergent. If not, the last step must be repeated with the starting hyperfine interactions parameters changed. After fitting is done, the software produces output text file with the data ready for importing in a graphical software like excell, origin, grapher etc..




Figure 16. The window of the MosK software with an experimental and calculated Mössbauer spectrum.

7 Results of measurements

The calibration of the spectrometer should be performed twice, i.e., before and after the measurement of the given sample. From these two measurements, the average value of the maximum velocity (in the last channel) is determined. This value is necessary for applying in the numerical procedure of the Mössbauer spectrum elaboration. The calibration spectrum (taken for α -iron) is presented in Figure 17.




Figure 17. Room-temperature calibration CEMS spectrum on α -iron foil.

From the numerical fit of the spectrum, the hyperfine interactions parameters are obtained. There are:

- isomer shift $IS = 0.01(1)$ mm/s (the value is relative to α -iron and should be zero);
- hyperfine magnetic field $B_{hf} = 33.06(2)$ T;
- quadrupole shift $QS = -0.003(5)$ mm/s (the value should be zero because iron has cubic lattice).

All the values are characteristic of α -iron.

The half width of the spectral lines $\Gamma = 0.16(2)$ mm/s is near the natural half width.

After calibration, the measurement for the proper sample is performed. In Figure 18, the CEMS spectrum for the steel sample Fe with 2.9 wt.% of Si is presented as an example.




Figure 18. Room-temperature CEMS spectrum of the sample Fe-2.9 wt.% Si with two components marked by red and green lines.

The main goal of the elaboration of the Mössbauer spectrum is to determine the number of the components and explain their origination. The fit is concerned as good when all the components are physically justified and the fit parameter χ^2 is near 1.



Funded by
the European Union

8 References

- [1] M. Przybylski, *From Three-to Two-Dimensional Magnetism of Fe*, habilitation dissertation, ed. by Wydawnictwa AGH, Kraków 1997.
- [2] M. Seberini, J. Sitek, *Surface Analysis of Iron-Containing samples by CEMS. An Overview of the Method*. Journal of electrical Engineering, 59 (5) 2008, 234-240.
- [3] B.V. Thosar, J.K. Srivastava, P.K. Iyengar, S.C. Bhargava, *Advances in Mössbauer Spectroscopy. Applications to Physics, Chemistry and Biology*. Elsevier Scientific Publishing Company, Amsterdam-Oxford-New York 1983.
- [4] Y. Garcia, J. Wang, T. Zhang, *Mössbauer Spectroscopy. Applications in Chemistry and Materials Science*. Wiley-VCH 2024.
- [5] I. Bibicu, *Mössbauer Spectroscopy Applied to Surfaces Study*. Rev. Roum. Chim. 57(9-10) 2012, 877-889.
- [6] W.R. Leo, *Techniques for Nuclear and Particle Physics Experiments*, 2nd ed. Springer-Verlag, Berlin, 1994.
- [7] M.F. Bent, B.I. Persson and D.G. Agresti, Computer Phys. Comm. 1, 1969, 67.
- [8] B.L. Christman and T.A. Tumolillo, Computer Phys. Comm. 2, 1971, 322.
- [9] W. Wilson and L.J. Swartzendruber, Computer Phys. Comm. 7, 1974, 151.
- [10] K. Kulcsar, D.L. Nagy and L. Pocs, Proc. Conf., Dresden, 1971, 594.
- [11] S.L. Ruby, *Mossbauer Effect Methodology*, vol. 8, ed. I.J. Gruverman and C.W. Seidel, Plenum Press, New York-London 1973.